Perfect™ Desensitizer

One drop hypersensitivity treatment

EN

For professional dental use only

Indications

- Perfect Desensitizer is used for treatment of dentin hypersensitivity.
- Prevents pain in exposed dentin areas not requiring restorations.
- Prevents hypersensitivity in dentin after preparation of teeth for following indirect restoration.

Properties

- Perfect Desensitizer is liquid desensitizer that quickly and effectively reduces pain associated with dentin hypersensitivity.
- Prevents intra-tubular fluid movement and reduces dental permeability by occlusion of the peripheral dentinal tubules. This inhibits the flow of fluid through the tubules which is the cause of sensitivity.
- The effect occurs immediately, effectively, simply, and fast with no need for mixing or curing.
- Perfect Desensitizer has been on the market since 2017 with no changes in formula. No complaints of any health-related or chemical performance issues were received.

Attention

Do not use in case of:

- Necessary precautions cannot be taken, or the specified application technique cannot be used
- Patient is allergic to any of the ingredients contained in Perfect Desensitizer.
- Pain is due to pulpitis.
- Pulp is exposed or the dentin adjacent to the pulp is perforated. Irritation of the pulp may be observed in the presence of deep cavities without suitable pulp protection. Those areas of cavities adjacent to the pulp should therefore be protected by a suitable cavity lining.

Caution!

- Perfect Desensitizer contains Glutaraldehyde, 2-Hydroxyethylmethacrylate.
- It is a local irritant and has toxic effects.
- Appropriate precautions should always be taken before using Perfect Desensitizer.
- Harmful by inhalation.
- Irritating to respiratory system.
- Irritating to skin.
- Risk of serious damage to eyes.
- Do not breathe vapors /spray.
- Avoid contact with skin.
- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- Wear suitable gloves.

- Before use, put on the protective goggles and cover the patient's eyes to protect against splashes of material.
- May cause sensitization by inhalation.
- May cause sensitization by skin contact.
- Perfect Desensitizer contains materials which: toxic / dangerous for the environment / toxic by inhalation and if swallowed / causes burns / may cause sensitization by inhalation and skin contact / very toxic to aquatic organisms / irritating to respiratory system and skin / may cause sensitization by skin contact. For complete symbols guide see below.

Composition

Essential Ingredients	% WT.
Water	<70%
Multifunctional monomers	<30%
Glutaraldehyde	<10%

Filler Content

The product does not contain filler.

Instruction for Use

Treatment of Hypersensitive Dentin:

The treatment of hypersensitive dentinal surfaces is only indicated in cases of severe pain in reactions to thermal, osmotic, or tactile stimuli. In cases of mild pain, alternative treatments such as desensitizing toothpastes should be used.

1. Clean the dentin (e.g., with cleansing pumice), under local anesthesia if necessary, and then rinse off with water.
2. Protect mucous membranes by using a rubber dam. Make sure that only Perfect Desensitizer comes into contact with the treated area.
3. Apply a small amount of Perfect Desensitizer to the dentinal surfaces using pellets or brushes and leave undisturbed for 30–60 seconds.
4. Dry the surface carefully by applying a gentle air flow until the fluid film has disappeared and the surface is no longer shiny.
5. Rinse thoroughly with water.
6. Tightly re-cap bottle after use.
7. If a single application of Perfect Desensitizer does not produce sufficient and sustained pain relief, the treatment can be repeated at the same or a subsequent appointment.

Treatment of Dentin Prior to Adhesive Bonding or Adhesive Cementation Procedures:

1. Protect mucous membranes by using a rubber dam.
2. Cover cavities and or areas near the pulp with suitable cavity linings before applying Perfect Desensitizer.
3. Remove moisture from the prepared dentin by applying a gentle air flow.
4. Apply a small amount of Perfect Desensitizer using pellets or brushes and leave undisturbed for 30–60 seconds.
5. Dry the dentinal surface carefully by applying a gentle air flow.
6. Rinse off thoroughly with water and apply suction.

7. Continue with Dentin Bonding or Cement of your choice and complete restoration.
8. Tightly re-cap bottle after use.

Storage and Disposal

- Do not store above 24°C (75°F).
- Keep product in original packaging.
- Keep away from direct sunlight.
- For optimum freshness, keep refrigerated.
- Do not freeze.
- To obtain SDS or IFU visit www.bjmlabs.com.
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.
- Empty product can be safely disposed after use in an intended regular waste container. The product does not require special disposal instructions or precautions and is harmless to the environment.

Symbols used on packaging

1. Consult instructions for use
2. Temperature limit
3. Date of manufacture
4. Use by date
5. Batch code
6. Medical device
7. Catalogue number
8. Irritant
9. Caution
10. Manufacturer

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B.J.M. Laboratories Ltd.
QUASER S.r.l.