
Dispensing / Distribution / Dispensación / Applicazione / Distribuição / Entnahme / Dozowanie / Distribuire / Παροχή

Application / Applicazione / Aplicación / Applicazione / Aplicação / Auftragen / Nakładanie / Aplicarea / Εφαρμογή

**Prima 1™**  
Light-cure self-etch 1 component dentin/enamel bonding agent

*For professional dental use only*

**EN INDICATIONS**

- Prima 1 is a self-etch one-component bonding agent for use in combination with adhesive restorations.
- Bonding of direct composite restorations.
- Used for Core Build-Up.
- Adhesive cementation of crown and bridges, inlays and onlays.

**PROPERTIES**

- Prima 1 is a 7th generation light-cure self-etch 1 component dentin and/or enamel bonding agent which combines Etch, Primer and Adhesive into a single component for ease of application.
- Prima 1 turns a multi-step procedure into a quick, easy solution that maximizes patient comfort while streamlining your working processes.
- Prima 1 forms an excellent bond with dentin and enamel and any self-cure or visible light-cure Composite, compomer restorative material and composite resin cements available on the market.
- Delivers high bonding values to dentin and enamel.
- Scanning electron micrographs indicates that Prima 1 application forms the essential hybrid zone with dentin.
- Prima 1 has been on the market since 2007 with no changes in formula. No complaints of any health-related or chemical performance issues were received.

**CAUTION!**

- Prima 1 contains polymerizable monomers, which may cause allergic contact dermatitis in susceptible individuals. Wash thoroughly with soap and water after contact. If skin sensitization occurs, or if known allergy to methacrylate resin exists, discontinue use.
- Polymerization may be triggered by ambient lighting or the dental operator light. Tightly recap the bottle after use.
- Prima 1 contains materials which: highly flammable / sensitive to light / irritating to eyes / irritating to respiratory system and skin / may cause sensitization by skin contact. For complete symbols guide see below.

**COMPOSITION**

Essential ingredients	% wt.	Essential ingredients	% wt.
Urethane dimethacrylate	<20%	Ethanol	<20%
Multifunctional monomers	<25%	Adhesive accelerator	<10%
Acetone	<40%	Photoinitiators	<5%

**FILLER CONTENT**  
The product does not contain filler.

**INSTRUCTIONS FOR USE**

**Tooth Preparation:**

- Prepare the tooth according to clinical requirements.
- Grind enamel before application. Prima 1 is less effective on uncut enamel.
- On sclerotic dentin, etch with Q-Etch or other 32% - 40% phosphoric acid etching gel available on the market for 30 seconds.
- Rinse with water and air dry.
- Isolate the cavity (the use of rubber dam is recommended).
- If indicated, place subbase.

**Note:** Do not use directly on exposed pulp tissue.

**Dispensing:**

- Shake bottle briefly prior to dispensing.
- Dispense Prima 1 into a dappen dish.
- Tightly re-cap bottle immediately after dispensing.
- Apply adhesive immediately after dispensing (within 3 minutes).

**Application of Prima 1:**

- Immediately after dispensing, soak the applicator tip or brush and apply ample amount to the entire cavity surface. Start with the enamel and move to the dentin. Apply 2 additional coats of Prima 1 slowly and consecutively.
- Leave undisturbed for 30 seconds. Slight agitation during the dwell time may improve demineralization and diffusion. Keep the coated surface free of contamination.
- Spread Prima 1 carefully with a gentle air flow to remove solvents, until no visible movement of the resin is noted. Continue air-drying for a few seconds with a stronger stream of air. **Note:** A strong initial air blast will scatter the bonding agent, resulting in poor adhesion.
- Surface should be visibly glossy, both after application of Prima 1 and after evaporation of the solvents. Ensure coverage of the entire cavity surface. If the cavity surface does not appear glossy, apply additional coats of Prima 1 and proceed as described above. Light-cure Prima 1 for 20 seconds with a very light-curing unit.
- Restore with light-cure or self-cure composite resin of your choice or cement restoration.

**STORAGE AND DISPOSAL**

- Do not store above 24°C (75°F).
- Keep product in original packaging.
- Keep away from direct sunlight.
- For optimum freshness, keep refrigerated.
- Do not freeze.
- To obtain SDS or IFU visit [www.bjmlabs.com](http://www.bjmlabs.com).
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.
- Empty product can be safely disposed after use in an intended regular waste container. The product does not require special disposal instructions or precautions and is harmless to the environment.

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