Prima 2000™

Single step bonding agent

For professional dental use only

EN INDICATIONS

- Direct bonding of composite restorations.
- Direct bonding of Core Build-Up.
- Adhesive cementation of crowns and bridges, inlays and onlays.

PROPERTIES

- 5th generation dentin and/or enamel bonding agent.
- Combines primer and adhesive into a single component for ease of application.
- Compatible with all common dental etching materials and forms an excellent bond with dentin and enamel.
- Compatible with any self-cure or light-cure composites, compomer restorative materials, and composite resin cements available on the market.
- Studies have demonstrated that Prima 2000 can deliver more than 35 MPa shear bond strength to etched enamel and more than 20 MPa shear bond strength to etched dentin.
- Scanning electron micrographs indicates that Prima 2000 application forms the essential hybrid zone with dentin.
- marketed since 1997 with no changes in formula. No complaints of any health-related or chemical performance issues were received.

CAUTION!

- Requires a calcium hydroxide liner/base when working near the pulp!
- Contains polymerizable monomers which may cause allergic contact dermatitis in susceptible persons. Wash thoroughly with soap and water after contact. If skin sensitization occurs, or if known allergy to methacrylate resin exists, discontinue use.
- Polymerization may be triggered by ambient lighting or the dental operator light. Tightly recap the bottle after use.
- Contains materials which: highly flammable / sensitive to light / irritating to eyes / irritating to respiratory system and skin / may cause sensitization by skin contact. For complete symbols guide see below.

COMPOSITION

Essential ingredients	% wt.	Essential ingredients	% wt.
Multifunctional methacrylates	<30%	Acetone	<50%
Multifunctional monomers	<10%	Ethanol	<10%
Urethane dimethacrylate	<10%	Adhesive accelerator	<5%
Photoinitiators	<1%	Filler	<1%



FILLER CONTENT

The filler system comprises:
0.5-1 weight-% or 0.1-0.5 volume-%.
The particle size of the fillers range: <6 um.

INSTRUCTIONS FOR USE

Tooth Preparation:

Prepare the tooth according to clinical requirements.

Etching / Mordançaço / Grabado / Mordenzatura / Condicionamento / Ätzen: Wytrawianie / Gravare / Αδρονολογή.

Application / Application / Aplicación / Applicazione / Aplicação / Auftragen / Nakładanie kleju / Aplicarea / Εφαρμογή του Prima 2000.

Etching:

1. Etch dentin for 5-10 seconds. For enamel, etch for 20-30 seconds using Q-Etch for both enamel and dentin.
2. Rinse with a spray of oil-free water for 5-10 seconds.
3. Blot dry excess moisture from the surface. Verify that the surface has no visible moisture.

Application of Prima 2000:

1. Dispense Prima 2000 directly onto a brush or a microbrush by holding the brush tip directly under the nozzle of the tilted bottle.
2. Apply ample amount of Prima 2000, leave undisturbed for 15 seconds.
3. Remove solvents with a gentle air flow for at least 15 seconds or until no visible movement of the resin is noted. If the surface appears partially coated, apply additional material and dry again.
4. Tightly re-cap bottle after use.
5. Light-cure Prima 2000 for 20 seconds with a light-cure unit.

If the restoration retention area lies primarily within dentin, apply a second layer of Prima 2000, remove solvents immediately and light cure for 20 seconds.

6. For light-inaccessible procedures, Prima 2000 can be self-cured by utilizing the Auto-Cure Activator. Mix equal number of drops of Prima 2000 and Auto-Cure Activator in a clean mixing well and proceed with step No.2.

7. Restore with light-cure or self-cure composite resin of your choice or cement restoration with composite resin cement.

STORAGE AND DISPOSAL

- Do not store above 24°C (75°F).
- Keep product in original packaging.
- Keep away from direct sunlight.
- For optimum freshness, keep refrigerated.
- Do not freeze.
- To obtain SDS or IFU visit www.bjmlabs.com.
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.
- Empty product can be safely disposed after use in an intended regular waste container. The product does not require special disposal instructions or precautions and is harmless to the environment.

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