Prima Quick™

Prima Quick Prime application
Application / Aplicación /
Aplicazione / Aplicação /
Auftragen / Nakładanie kleju /
Aplicarea / Εφαρμογή tou.

Prima Quick™

Universal, no-rinse, self-priming etchant adhesive system

For professional dental use only

EN INDICATIONS

- Direct composite restorations.
- Adhesive cementation of crowns and bridges, inlays and onlays.

PROPERTIES

- Prima Quick is a 6th generation self-etch bonding system, which delivers high bonding values to both dentin and enamel.
- Prima Quick consists of an acidic water based self-etch Primer and a light-cure Adhesive.
- It is used to bond composite resin filling materials, to tooth structure, any self-cure or light-cure composite, compomer restorative materials and composite resin cements available on the market.
- Prima Quick is utilizing a moderately acidic Primer. Unlike total-etch, total-bond techniques, the self-etch Primer does not open the dentinal tubules completely. The smear layer is solubilized, and due to the high hydrophilic properties of the Primer, it is then able to penetrate to the tubules and peritubular dentin, forming resin tags. In the case of enamel, the Primer produces a significant pattern with enhanced surface area, leading to improved enamel bonding.
- Prima Quick has been on the market since 2008 with no changes in formula. No complaints of any health-related or chemical performance issues were received.

CAUTION!

- Prima Quick contains polymerizable monomers which may cause allergic contact dermatitis in susceptible individuals. Wash thoroughly with soap and water after contact. If skin sensitization occurs, or if you know allergy to methacrylate resin exists, discontinue use.
- Polymerization may be triggered by ambient lighting or the dental operatory light. Tightly recap the bottle after use.
- Prima Quick contains materials which: sensitive to light/ highly flammable / irritating to eyes / harmful if swallowed / irritating to respiratory system and skin / may cause sensitization by skin contact / repeated exposure may cause skin dryness or cracking / vapors may cause drowsiness and dizziness. For complete symbols guide see below.

COMPOSITION

Essential ingredients – PRIME	% wt.
Multifunctional monomers	<20%
Methacrylated phosphoric acid esters	<5%
Adhesive accelerator	<10%
Water	<40%
Ethanol	<30%
Acetone	<10%

FILLER CONTENT

The filler system comprises:
0.5-1 weight-% or 0.1-0.5 volume-%.
The particle size of the fillers range: <6 µm.

Essential ingredients – BOND	% wt.
Multifunctional methacrylates	<30%
Multifunctional monomers	<10%
Urethane dimethacrylate	<10%
Photoinitiators	<1%
Acetone	<50%
Ethanol	<10%
Adhesive accelerator	<5%
Filler	<1%

Prima Quick Bond application
Application / Aplicación /
Aplicazione / Aplicação /
Auftragen / Nakładanie kleju /
Aplicarea / Εφαρμογή tou.

INSTRUCTIONS FOR USE

Application of Prima Quick Prime:

1. Apply ample amount of Prima Quick Prime on the prepared, air-dried or blotted tooth surface using a sponge or disposable brush tip. Leave undisturbed for 15 seconds.
2. Remove solvents with a blast of air for at least 15 seconds.
3. Tightly re-cap bottle after use.

Application of Prima Quick Bond:

1. Dispense Prima Quick Bond directly onto a brush or a microbrush by holding the brush tip directly under the nozzle of the tilted bottle.
2. Apply ample amount of Prima Quick Bond, leave undisturbed for 15 seconds.
3. Remove solvents with a gentle air flow for at least 15 seconds or until no visible movement of the resin is noted. If the surface appears partially coated, apply additional material and dry again.
4. Tightly re-cap bottle after use.
5. Light-cure Prima Quick Bond for 20 seconds with a light-cure unit.
If the restoration retention area lies primarily within dentin, apply a second layer of Prima Quick Bond, remove solvents immediately and light cure for 20 seconds.
6. For light-inaccessible procedures, Prima Quick Bond can be self-cured by utilizing the Auto-Cure Activator. Mix equal number of drops of Prima Quick Bond and Auto-Cure Activator in a clean mixing well and proceed with step No.2.
7. Restore with light-cure or self-cure composite resin of your choice or cement restoration with composite resin cement.

STORAGE AND DISPOSAL

- Do not store above 24°C (75°F).
- Keep product in original packaging.
- Keep away from direct sunlight.
- For optimum freshness, keep refrigerated.
- Do not freeze.
- To obtain SDS or IFU visit www.bjmlabs.com.
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.
- Empty product can be safely disposed after use in an intended regular waste container. The product does not require special disposal instructions or precautions and is harmless to the environment.

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