Q-Core™ / Q-Core™ Syringable

Dual-cure composite core build-up material

For professional dental use only

EN INDICATIONS

- Q-Core is a dual-cure, radiopaque resin composite material designed for the fabrication of core build-ups.
- Q-Core may be used with any Bis-GMA compatible bonding agent.
- Used for core build-ups in vital and non-vital teeth.

PROPERTIES

- Q-Core is a restorative material.
- Q-Core classification according to ISO 4049 is Type 2, Class 3.
- Easy to use.
- Provides a consistent, homogeneous mix.
- Q-Core is a dual-cure composite – it has a self-curing mechanism which can also be light cured. The self-curing assures that the material can be quickly bulk filled.
- Dense void-free core: cuts like dentin.
- Q-Core has a radiopacity value above 4.0 mm Aluminium. (Radiopacity is defined as the ability of a material to be visible in x-ray photographs).
- Superior compressive strength for long-lasting restorations.
- Q-Core is available in 3 shades: A3, White and Blue. To see shade guide visit www.bjmlabs.com.
- Q-Core has been on the market since 2004 with no changes in formula. No complaints of any health-related or chemical performance issues were received.

CAUTION!

- Q-Core contains polymerizable monomers which may cause allergic contact dermatitis in susceptible patients. Wash thoroughly with soap and water after contact. If skin sensitization occurs, or if known allergy to methacrylate resin exists, discontinue use.
- Q-Core contains materials which: sensitive to light / irritating to eyes / irritating to respiratory system and skin / may cause sensitization by skin contact. For complete symbols guide see below.

COMPOSITION

Essential ingredients – BASE	% wt.	Essential ingredients – CATALYST	% wt.
Triethylenglycoldimethacrylate	<20%	Triethylenglycoldimethacrylate	<20%
Bisphenol A-glycidyl methacrylate	<20%	Bisphenol A-glycidyl methacrylate	<20%
Polymerization activators	<5%	Benzoyl Peroxide	<1%
Photoinitiators	<1%	Glass filler	<60%
Glass filler	<60%	Fumed Silica	<5%
Fumed silica	<5%		
Pigments	<1%		

FILLER CONTENT

The filler system comprises glass fillers and modified silica: 60-65 weight-% or 15-25 volume-%.

The particle size of the fillers ranges: <6 um.

INSTRUCTIONS FOR USE

- Note:**
- Do not use in presence of eugenol-based materials; complete setting of Q-Core may not take place.**

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- A calcium hydroxide or glass ionomer liner may be used to isolate Q-Core from eugenol sources.**
- If refrigerated, allow to return to room temperature before use. Refrigeration will increase viscosity, making application more difficult, and will retard setting time.**
- Tips are for single patient use only to prevent cross-contamination between patients.**

- Prepare the tooth according to clinical requirements.
- Isolate the tooth from oral fluids. A rubber dam is recommended.
- Etch the entire preparation site with Q-Etch or other 32%-40% phosphoric acid etching gel available on the market for 15 seconds or apply Prima Quick Prime and leave undisturbed for 15 seconds. When using another manufacturer's etchant, follow manufacturer instructions.
- Rinse off etchant with a water spray for at least 15 seconds. Etched surface should appear frosted. When using Prima Quick Prime, remove solvents with a blast of air for at least 10 seconds.
- Mix equal drops of Prima 2000 or Prima Quick Bond with the Auto-Cure Activator and place on the prepared tooth surface. When using another manufacturer's bonding agent, follow instructions for use with dual-cure composites.
- Allow mixture to sit for 15 seconds and then remove residual solvents with a gentle blast of oil-free air. Repeat application of Prima 2000 / Prima Quick Bond mixture at least once on the prepared cavity tooth surface. It is not necessary to wait 15 seconds between each additional application; however, solvent evaporation is necessary between applications.
- Light-cure for 20 seconds after final application of bonding agent is placed.
- Bleed cartridge/syringe before installation of the automix dispensing tip/automix syringe mix tip. Mount a new, automix dispensing tip onto the cartridge/automix syringe mix tip onto the dual-barrel syringe and bleed the cartridge again/syringe. This ensures complete flow of Base and Catalyst components.
- Place into base of prepared tooth an adequate amount of Q-Core to complete build-up.
- Do not disturb Q-Core during the gelation stage. Light-cure for at least 40 seconds. Check for complete hardness.
- Tightly re-cap syringe or cartridge after use.
- Time indications:
 - The working and setting times depend on the ambient and oral temperature.
 - The times shown are based on conditions relevant for practice.

	min:sec
Working time at ambient temperature	>01:30
Intra-oral setting time	<05:00
Light curing in a light-wave range of 390 – 470 nm <ul style="list-style-type: none">Compatible with all curing units available on the market	00:40

STORAGE AND DISPOSAL

- Do not store above 24°C (75°F).
- Keep product in original packaging.
- Keep away from direct sunlight.
- For optimum freshness, keep refrigerated.
- Do not freeze.
- If refrigerated, allow the syringe to reach room temperature.
- To obtain SDS or IFU visit www.bjmlabs.com.
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.
- Empty product can be safely disposed after use in an intended regular waste container. The product does not require special disposal instructions or precautions and is harmless to the environment.

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