Q-Glass™

Glass-ionomer self-cure luting cement

For professional dental use only

EN INDICATIONS

- Q-Glass is a glass ionomer cement for permanent cementation of crowns and bridges.
- Definitive sealing of inlays, crowns, and bridges.

PROPERTIES

- Easy mix.
- High strength.
- Snap setting.
- Easy removal of excess material.
- Can also be used as a liner and base.
- Q-Glass has been on the market since 2000 with no changes in formula. No complaints of any health-related or chemical performance issues were received.

CAUTION!

- Q-Glass Liquid and mixed Q-Glass cement paste are corrosive. Avoid contact with skin, oral soft tissues, and eyes. In case of contact with oral tissue or skin, remove immediately with alcohol-soaked cotton wool and rinse thoroughly with water. In case of contact with eyes, immediately rinse thoroughly with water and seek medical attention.
- Do not mix Q-Glass powder or liquid with any other glass ionomer-based component system. Direct contact with eugenol-based products is to be avoided. Eugenol inhibits the setting of glass ionomer sealers.
- Q-Glass contains materials which: irritating to eyes / irritating to respiratory system and skin. In the event of an allergic reaction, consult a specialist. For complete symbols guide see below.

COMPOSITION

Essential ingredients – POWDER	% wt.	Essential ingredients – LIQUID	% wt.
Strontium aluminosilicate glass powder	<90%	Polyacrylic acid	<70%
Polyacrylic acid	<10%	Tartaric acid	<15%
Hydroxyapatite	<10%	Water	<25%

FILLER CONTENT IN MIXED PASTE

The filler system comprises glass fillers and modified silica:

80-90 weight-% or 30-40 volume-%.

The particle size of the fillers range: <6 µm.

INSTRUCTIONS FOR USE

- Prepare the tooth according to the usual protocol.

- Clean with water and pumice. Rinse and dry without draining so as not to alter the “smear layer”. For pulp capping, use a calcium hydroxide cement.
- Clean and dry the intrados of the prosthetic part.
- Before dispensing, fluff the powder bottle by shaking it up and down for 5 seconds and shake the liquid bottle. For accurate dispensing of liquid drops, hold the bottle vertically upside down. Allow air bubbles to escape before dripping. If the nozzle is clogged, clear with damp gauze. Discard drops with air bubbles.
- Mix the necessary quantity of cement. The recommended dosage is 1 spoon of powder for 3 drops of liquid, which corresponds to 0.16 gr powder and 0.12 gr liquid.
Powder/Liquid ratio (g/g): 0,16 / 0,12
- Apply one level scoop of powder and 3 drops of liquid on a mixing pad.
- Using a clean spatula introduce the powder into the liquid and mix briskly for 20 seconds. If larger quantities are required, divide the powder into 2 equal portions:
Mix one portion with the liquid for 5 seconds. Add the second portion of powder and mix for 15 seconds.
- Apply a 1 mm layer of cement to the inside surface of the prosthetic and position it within 30 seconds after mixing has ceased.
Note: High temperatures short working time.
- Setting time is 5 minutes from the beginning of mixing. Once it has hardened, remove the excess material, and protect the prosthetic joint with a polish or unfilled polymerizable resin.
- Tightly re-cap the bottles after use.

Time indications:

The working and setting times depend on the ambient and oral temperature.

The times shown are based on conditions relevant for practice.

	min:sec
Mixing time	00:20
Working time	2:10-2:30
Setting time (at 37°C)	4:00-4:30

STORAGE AND DISPOSAL

- Do not store above 24°C (75°F).
- Keep product in original packaging.
- Keep away from direct sunlight.
- For optimum freshness, keep refrigerated.
- Do not freeze.
- If refrigerated, allow the syringe to reach room temperature.
- To obtain SDS or IFU visit www.bjmlabs.com.
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.
- Empty product can be safely disposed after use in an intended regular waste container. The product does not require special disposal instructions or precautions and is harmless to the environment.

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