High-Q-Bond SE™

High-Q-Bond SE™

Self-etch, self-adhesive dual-cure permanent resin cement

For professional dental use only.

Patient Population: Any human undergoing dental treatment related to the procedures listed below under "Indications". The definitive patient group is diagnosed on a case-by-case basis by a licensed dental professional and there are no exclusions for specific patient groups.

EN INDICATIONS

- Permanent cementation of crowns and bridges, inlays and onlays, posts and cores, ceramic crowns and Maryland bridges.
- Metall alloys adhesion opaquer.
- Porcelain repair.
- Proven Color Stability over time.

PROPRIETIES

- High-Q-Bond SE is a self-etch, self-adhesive dual-cure, radiopaque permanent resin cement and does not require etching, priming or bonding.
- High-Q-Bond SE classification according to ISO 4049 is Type 2, Class 3.
- High-Q-Bond SE has a radiopacity value above 2.5 mm Aluminium. Radiopacity is defined as the ability of a material to be visible on x-ray photographs.
- High-Q-Bond SE forms ideal bond with metal alloys, porcelain, ceramic, dentin, enamel, amalgam.
- High-Q-Bond SE has a chemical affinity for metal.
- High-Q-Bond SE is a dual-cure cement – it has a self curing mechanism which can also be light cured for an immediate set and extra strength at the margins.
- Light curing High-Q-Bond SE gives restorations immediate stability and easy clean up of excess cement.
- High-Q-Bond SE features high retention and total margin integrity.
- High-Q-Bond SE Automix syringe saves application time and guarantees consistent mix.
- High-Q-Bond SE is available in 3 shades: A2, White and Translucent. To see shade guide visit www.bjmlabs.com.
- High-Q-Bond SE has been on the market since 2008 with no changes in formula. No complaints of any health-related or chemical performance issues were received.

CONTRAINDICATIONS

- Not to be used on exposed pulp.
- The use of the medical device is contraindicated if the patient is known to be allergic or sensitive to any of its ingredients.
- Direct contact with eugenol-based products is to be avoided. Eugenol inhibits the setting of resin-based cements.

CAUTION!

- High-Q-Bond SE contains polymerizable monomers which may cause allergic contact dermatitis in susceptible patients and damage the pulp. Avoid contact with skin, eyes and soft tissue. Wash thoroughly with water after contact. If skin sensitization occurs or if known allergy to methacrylate resin exists, discontinue use.
- High-Q-Bond SE contains materials which: sensitive to light / irritating to eyes / harmful if swallowed / irritating to respiratory system and skin / may cause sensitization by skin contact. For complete symbols guide see below.

COMPOSITION

Essential ingredients – BASE	% wt.	Essential ingredients – CATALYST	% wt.
Bisphenol A-glycidyl methacrylate	<40%	Bisphenol A-glycidyl methacrylate	<40%
Triethyleneglycoldimethacrylate	<20%	Triethyleneglycoldimethacrylate	<20%
Polymerization activators	<1%	2-Hydroxyethylmethacrylate	<10%
Photoinitiators	<1%	4-Methacryloxyethyltrimellitic anhydride	<5%
2-Hydroxyethylmethacrylate	<10%	Benzoyl Peroxide	<1%
Glass filler	<40%	Glass filler	<50%
Fumed silica	<5%	Fumed silica	<5%

FILLER CONTENT

The filler system comprises glass fillers and modified silica: 40-50 weight-% or 10-20 volume-%.

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INSTRUCCIONES FOR USE

Nota:

- Tips are for single patient use only to prevent cross-contamination between patients.
- Use one Detachable Manual Adapter per one syringe lifetime.
- Remove the provisional restoration. **Note:** Use only temporary cements that are eugenol free. The use of eugenol containing temporary materials can inhibit the polymerization process of the High-Q-Bond SE.
- Rinse thoroughly.
- Air dry the preparation (do not desiccate).
- Try in the restoration to ensure proper fit.
- Automix (When using a dual-barrel syringe with automix syringe mix tip):** Bleed the syringe before installation of the automix syringe mix tip. Mount a new, automix syringe mix tip onto the dual-barrel syringe and bleed the syringe again. This ensures complete flow of Base and Catalyst components.
- High-Q-Bond SE Automix Syringe with Detachable Manual Adapter:** a) Bleed the syringe before installation of the Detachable Manual Adapter. Mount a new manual adaptor onto the dual-barrel syringe and bleed the syringe again. This ensures complete flow of Base and Catalyst components. b) Using a clean spatula or mixing stick, mix the two components into a homogenous paste.
- Dispensing the material can be done as follows:
 - For Inlay / Onlay** - Mount the syringe intra oral tip onto the mixing tip and dispense the cement directly, covering all surfaces. Seal the restoration gently into the preparation allowing the cement to flow from all sides.
 - For Crown** - Dispense the cement directly into the restoration, covering all surfaces. Seal the restoration gently onto the preparation allowing the cement to flow from all margins.
 - For Post** - Mount the supplied XFX intra oral tip onto the mixing tip. Syringe the cement into the canal. Seal the post and vibrate slightly avoiding possibility of trapped air. Allow the cement to slowly flow from the canal space.
- When the restoration is properly seated, remove all excess cement. Excess cement is best removed in its gel state, after reaching the material working time, with a scaler or explorer. After removal of excess cement, light cure all surfaces including margins for 20 seconds.
- For the material time indications, see table below.
- Time indications:
 - The working and setting times depend on the ambient and oral temperature.
 - The times shown are placed on conditions relevant for practice.

	min:sec
Working time at ambient temperature	>01:30
Intra-oral setting time	<05:00
Light curing in a light-wave range of 390 – 470 nm Compatible with all curing units available on the market	00:20

- Tightly re-cap syringe after use.
- Use finishing strips for proximal areas. Polish margins with polishing discs, cups, or points.
11. Check occlusion.

Mixing options:

Note: working and setting times complies to the mentioned above range without dependence on mixing options.

Mixing options	Ratio
Automix - 1 Automix syringe	1:1 - Equal amount of Base and Catalyst
Hand Mix - 1 Automix syringe with Detachable Manual Adapter	1:1 - Equal amount of Base and Catalyst

STORAGE AND DISPOSAL

- Do not store at temperatures above 25°C (75°F).
- Keep product in original packaging.
- Keep away from direct sunlight.
- For optimum freshness, keep refrigerated.
- Do not freeze.
- If refrigerated, allow the syringe to reach room temperature.
- To obtain SDS or IFU visit www.bjmlabs.com.
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.
- Empty product can be safely disposed after use in an intended regular waste container. The product does not require special disposal instructions or precautions and is harmless to the environment.

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