## Zirconite™

### Dual-cure permanent adhesive resin cement designed for zirconia

For professional dental use only.

**Patient Population:** Any human undergoing dental treatment related to the procedures listed below under "Indications". The definitive patient group is diagnosed on a case-by-case basis by a licensed dental professional and are there are exclusions for specific patient groups.

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- Permanent Cementation of Ceramic and Zirconia restorations.
- Proven Color Stability over time.

#### PROPERTIES

- A novel dual-cure permanent adhesive resin cement specially formulated for Zirconia cementation.
- ISO 4049 classification: Type 2, Classe 3.
- Self-etch, self-adhering, dual cured, radiopaque permanent resin cement and does not require etching, priming, or bonding.
- Radiopacity value above 2.5 mm Aluminium. (Radiopacity is defined as the ability of a material to be visible in x-ray radiographs.)
- Chemical affinity for metal and Zirconia.
- Dual-cure cement – it has a self curing mechanism which can also be light cured for an immediate set and extra strength at the margins.
- Features high retention and total margin integrity.
- Administration via Automix syringe allows application time and guarantees consistent mix.
- Available in 3 shades: Dentin, White, Translucent. To see shade guide visit [www.bjmlabs.com](http://www.bjmlabs.com).
- Zirconite has been on the market since 2009 with no changes in formula. No complaints of any health-related or chemical performance issues were received.

#### CONTRAINDICATIONS

- Not to be used on exposed pulp.
- The use of the medical device is contraindicated if the patient is known to be allergic or sensitive to any of its ingredients.
- Direct contact with eugenol-based products is to be avoided. Eugenol inhibits the setting of resin-based cements.

#### CAUTION!

- Contains polymerizable monomers which may cause allergic contact dermatitis in susceptible patients and damage the pulp. Avoid contact with skin, eyes, and soft tissue. Wash thoroughly with water after contact.
- If skin sensitization occurs or if known allergy to methacrylate resin exists, discontinue use.
- Contains materials which are sensitive to light / irritating to eyes / irritating to respiratory system and skin / may cause sensitization by skin contact. For complete symbols guide see below.
- Any serious incident that occurs in relation to this device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established. You may also contact our EU Authorised Representative: Obelis S.A. (contact details below).

#### COMPOSITION

Essential ingredients - BASE	% wt.	Essential ingredients - CATALYST	% wt.
Composition of acrylic resins	<20%	Composition of Acrylic Resins	<30%
2-Hydroxyethylmethacrylate	<10%	Triethylenglycoldimethacrylate	<20%
Triethylenglycoldimethacrylate	<10%	Sodium p-toluensulfinate	<10%
Photoinitiators	<5%	Fumed Silica	<5%
Copper (II) chloride dihydrate (3% aqueous solution)	<1%	Glass filler	<60%
10-methacryloyloxydecyl dihydrogen phosphanate	<10%		
Fumed silica	<5%		
Glass filler	<60%		

#### FILLER CONTENT

The filler system comprises glass fillers and modified silica:

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#### INSTRUCTIONS FOR USE

##### Notes:

- This are for single patient use only to prevent cross-contamination between patients.
- Use one Detachable Manual Adaptor per one syringe lifetime.
- 1. Remove the provisional restoration.
- Note: Use only temporary cements that are eugenol free. The use of eugenol-containing temporary luting materials can inhibit the polymerization process of Zirconite.
- 2. Rinse thoroughly.
- 3. Air dry the preparation (do not desiccate).
- 4. Try in the restoration to ensure proper fit.
- 5. Automix (When using a dual-barrel syringe with Automix syringe mix tip): Bleed the syringe before installation of the automix syringe mix tip. Mount a new, automix syringe mix tip onto the dual-barrel syringe and bleed the syringe again. This ensures complete flow of Base and Catalyst components.
- Hand Mix (When using a dual-barrel syringe with Detachable Manual Adapter): a) Bleed the syringe before installation of the Manual Adapter. Mount a new manual adaptor onto the dual-barrel syringe and bleed the syringe again. This ensures complete flow of Base and Catalyst components. b) Using a clean spatula or mixing stick, mix the two components into a homogenous paste.
- 6. Dispensing the material can be done as follows:
  - For Inlay/Onlay - Mount the expander into oral tip onto the mixing tip and dispense the cement directly into the cavity preparation covering all surfaces. Seat the restoration gently onto the preparation allowing the cement to flow from all sides.
  - For Crown - Dispense the cement directly into the restoration or tooth preparation covering all surfaces. Seat the restoration gently onto the preparation allowing the cement to flow from all margins.
- Note: For the material time indications, see table below.
- 7. When the restoration is properly seated, remove all excess cement. Excess cement is best removed in its gel state with a scaler or explorer.
- 8. Tightly re-cap syringe after use.
- 9. Time indications:

The working and setting times depend on the ambient and oral temperature.

The times shown are based on conditions relevant for practice.

	min:sec
Working time at ambient temperature	>01:30
Intra-oral setting time	<05:00
Light curing in a light-wave range of 390 – 470 nm	00:20
Compatible with all curing units available on the market	

- 10. After removal of excess cement, light cure all surfaces including margins for 20 seconds.
- 11. Use of air strips for proximal areas. Polish margins with polishing discs, cups, or points.
- 12. Check occlusion.

#### Mixing options:

Note: working and setting times complies to the mentioned above range without dependence on mixing options.

Mixing options	Ratio
Automix - 1 Automix syringe	1:1 - Equal amount of Base and Catalyst
Hand Mix - 1 Automix syringe with Detachable Manual Adapter	1:1 - Equal amount of Base and Catalyst

#### STORAGE AND DISPOSAL

- Do not store above 24°C (75°F).
- Keep product in original packaging.
- Keep away from direct sunlight.
- For optimum freshness, keep refrigerated.
- Do not freeze.
- If refrigerated, allow the syringe to reach room temperature.
- To obtain SDS or IFU visit [www.bjmlabs.com](http://www.bjmlabs.com).
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.
- Empty product can be safely disposed after use in an intended regular waste container. The product does not require special disposal instructions or precautions and is harmless to the environment.

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