# Q TEMP™

## EN NON-EUGENOL, RESIN BASED TEMPORARY CEMENT CONTAINING FLUORIDE, CHLORHEXIDINE AND POTASSIUM NITRATE

### INDICATIONS FOR USE

Q-Temp Temporary Cement is recommended for the cementation of temporary crowns, bridges, inlays and onlays.

### CONTRAINDICATIONS

Q-Temp Temporary Cement is not intended as a temporary filling material nor is it to be used for pulp capping.

### INSTRUCTIONS FOR USE

1. Complete tooth preparation and fabrication of the temporary restoration before extruding the mixed cement.
2. Dry and isolate the prepared tooth. Be careful not to desiccate or over-dry the preparation.
3. Remove cap from double barrel automix syringe and set aside. Extrude a very small portion of the unmixed material to insure unobstructed flow from both the base and catalyst components of the respective barrels.
4. Attach automix dispensing tip and lock in place with a quarter turn.

5. Extrude a thin layer of the automixed cement on to the cleaned, inside surface of the temporary restoration. You have approximately 45 seconds of working time after extruding the mixed cement.
6. Firmly seat the temporary restoration in place. Initial gel-set occurs within 2 to 2.5 minutes after seating. Cleaning of excess cement should be initiated immediately after this gel-set phase. 7) Final cement set occurs between 4 and 5 minutes after seating of the temporary restoration.
7. Remove used mixing tip from double barrel automix syringe and replace with original cap. If cap is misplaced, the used tip can be left in place until the next use.

### HAND-MIXING DIRECTIONS

1. Complete tooth preparation and fabrication of the temporary restoration before extruding the mixed cement.
2. Dry and isolate the prepared tooth. Be careful not to desiccate or over-dry the preparation.
3. Extrude the appropriate amount of base and catalyst onto the selected mixing surface.

4. Using a clean spatula or mixing stick, mix the two components into a homogenous paste.
5. Using the mixing stick or brush, apply a thin layer of the mixed cement on to the cleaned, inside surface of the restoration. You have approximately 45 seconds of working time from start of mix.
6. Firmly seat the temporary restoration in place. Initial gel-set occurs within 2 to 2.5 minutes after seating. Cleaning of excess cement should be initiated immediately after this gel-set phase.
7. Final cement set occurs between 4 and 5 minutes after seating of the temporary restoration.

### PRECAUTIONS

Though rare, some individuals may be sensitive to this product. If an allergic reaction occurs, discontinue use, remove from patient and refer patient to a physician. Important Application Information: Only a small bead of Q-Temp around the margin is necessary to securely hold your provisional. This is especially true with highly retentive temporaries. Using less cement saves money and allows for easy retrievals!

### ESSENTIAL INGREDIENTS

#### Base & Catalyst:

- Triethyleneglycoldimethacrylate.
- Urethane dimethacrylate.
- Polymerization activators.
- Fumed Silica.

#### STORAGE

- Do not store above 24°C (75°F).
- For optimum freshness, keep refrigerated.
- Do not freeze.
- If refrigerated, allow the syringe to reach room temperature.
- To obtain SDS visit [www.bjmlabs.com](http://www.bjmlabs.com).
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.

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