Prima 1™

EN LIGHT CURED SELF ETCHING 1 COMPONENT DENTIN/ENAMEL BONDING AGENT

PROPERTIES

Prima 1 is a light cured self etching 1 component dentin/enamel bonding agent which combines etch, primer and adhesive into a single component for ease of application. Now you can save time and efforts by etching, priming and bonding in a single step. No mixing. No mess. No worries.

Prima 1 turns a multi-step procedure into a quick, easy solution that maximizes patient comfort while streamlining your working processes.

Prima 1 utilizes the latest advancements in nanotechnology and dendritic polymer incorporated chemistry and forms an excellent bond with dentin and enamel, any self-cure or visible light-cured Composite, Compomer restorative material and with Composite Resin Cements.

Scanning electron micrographs indicate that Prima 1 application forms the essential hybrid zone on dentin. Now you can be confident that your bonding process will be perfect every time.

INDICATIONS

- Direct composite restorations
- Core build up
- Adhesive cementation of crown & bridge, including inlays and onlays

CAUTION

Prima 1 contains polymerizable monomers which may cause allergic contact dermatitis in susceptible individuals. Wash thoroughly with soap and water after contact. If skin sensitization occurs, or if known allergy to methacrylate resin exists, discontinue use.

INSTRUCTIONS FOR USE

TOOTH PREPARATION

- Prepare the tooth according to clinical requirements.
- Rinse with water and air dry.
- Isolate the cavity (the use of rubber dam is recommended)
- If indicated, place subbase.

Do not use directly on exposed pulp tissue.

DISPENSING

Bottle application:

- Shake bottle briefly prior to dispensing.
- Dispense Prima 1 into a mixing well.
- Replace bottle cap immediately after dispensing. Apply adhesive immediately from dispensing. (Within 3 minutes)

APPLICATION OF PRIMA 1

- Immediately after dispensing, soak the applicator tip or brush and apply copious amount to the

entire cavity surface. Start with the enamel and move to the dentin. Apply 2 additional coats of Prima 1 slowly and consecutively. Cut/grind enamel before application. Prima 1 is less effective on uncut enamel. On sclerotic Dentin, etch with phosphoric acid for 30 seconds.

- Leave undisturbed for 30 seconds. Slight agitation during the dwell time may improve demineralization and diffusion. Keep the coated surface free of contamination.
- Spread Prima 1 carefully with a gentle stream of air to volatilize solvent, until no movement of the adhesive film is noticed. Continue air-drying for a few seconds with a stronger stream of air. The objective is to volatilize the solvent and water without removing the active ingredients from the surface. A strong initial airblast will scatter the bonding agent, resulting in poor adhesion.
- Surface should be visibly glossy, both after application of Prima 1 and after evaporation of the solvent. Ensure coverage of the entire cavity surface. If the cavity surface does not appear shiny, apply additional coats of Prima 1 and proceed as described above.
- Cure Prima 1 for 20 seconds with a halogen light-cure unit. Curing time assumes comparable intensity while use minimum 400 mW/cm². LED units require the same polymerization time. While use high intensity light curing units with an intensity > 1200

mW/cm², the curing time can be reduced by 25%. Low intensity light causes poor adhesion. Curing lights should be checked for proper output using a reliable light tester at periodic intervals. The emitting tip should be placed as close as possible to the resin surface.

- Restore with visible light cured or self-cured composite resin of your choice or cement restoration with Composite Resin Cement.

ESSENTIAL INGREDIENTS

- Acetone.
- Ethanol.
- Multifunctional dimethacrylates.
- Triéthylène glycol diméthacrylate.
- Multifunctional monomers.
- Adhesive accelerator.
- Photoinitiators.

STORAGE

- Do not store above 24°C (75°F).
- Keep away from direct sunlight.
- For optimum freshness, keep refrigerated.
- Do not freeze.
- Tightly recap after use.
- To obtain SDS visit www.bjmlabs.com.
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.

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