



"Improving Patient Care Through Research & Education"

Home > Clinical Evaluations > Evaluations

Premier Implant Cement

++++1/2



Premier Dental Products Co
888.670.6100 - 610.239.6000
www.premusa.com

Featured in:
January-February
2005
Volume 22, No. 01

Consultants' Comments

- ▶ "Very easy to mix and use."
- ▶ "Consistency is smooth and creamy."
- ▶ "Very retentive material."
- ▶ "I would like to see a selection of shades for use under anterior provisionals."
- ▶ "Include more syringes in the kit."
- ▶ "Include more mixing tips – the supply of tips is not adequate for the amount of material in the syringe."

Subscriber Toolbox

[READ FULL ISSUE](#)
[ADD TO MY ARCHIVE](#)
[PRINT THIS PAGE](#)

Description

Premier Implant Cement is a eugenol-free, elastic polymer cement indicated for use with long-term temporaries and implant-retained crowns that may require future removal. The cement has a two-stage cure – an initial gel-phase in 2.5 minutes for removal of excess cement and a rigid final set. It is a neutral, light pink gingival color. The kit includes a 5-ml dual-barrel, automix syringe, 10 mixing tips and one implant scaler instrument. Thirty consultants evaluated **Premier Implant Cement** in the cementation of 175 crowns.

Product Features

Premier Implant Cement was highly rated for its easy removal and clean-up of excess cement. The automix feature provides neat, quick mixing and easy placement into the restoration. The light viscosity of the cement allows complete seating of restorations without excessive pressure. Consultants favored the retentive quality of the material, although a few found it difficult to remove restorations when necessary. One hundred percent of the consultants found the cement to be equivalent to or better than their current cement. Sixty percent would switch, while ninetythree percent would recommend **Premier Implant Cement**.

Clinical Tip

- Incorporate a water-based lubricant into the mix to facilitate future removal of the crown.

More About: [Premier Dental Products Co, Cement - Temporary, Cement,](#)

CE0340

[FAQ](#) [Copyright Policy](#) [Privacy Policy](#) [Terms & Conditions](#) [Contact](#)

Dental Consultants, Inc. 3110 West Liberty, Ann Arbor, Michigan 48103
(800) 347-1330 (734) 665-2020 (734) 665-1648 fax info@dentaladvisor.com
© 2007 Dental Consultants, Inc.