



"Improving Patient Care Through Research & Education"

Home > Clinical Evaluations > Evaluations

Premier Implant Cement (abridged) +++++1/2

2007
PREFERRED
PRODUCT

Premier Dental Products Co
888.670.6100 - 610.239.6000
www.premusa.com

Consultants' Comments

- ▶ "Very easy to mix and use. Consistent smooth and creamy."
- ▶ "Very retentive material, but I would like to see a selection of shades for use on anterior provisionals."
- ▶ "Include more syringes in the kit."
- ▶ "Include more mixing tips – the supply of tips is not adequate for the amount of material in the syringe."

Preferred Products are selected by Dr. John Farah and the staff of THE DENTAL ADVISOR. These products are featured each year in the January-February issue of THE DENTAL ADVISOR.

Featured in:
January-February 2007
Volume 24,
No. 01

Subscriber
box

READ FULL
ISSUE

ADD TO MY
ARCHIVE

PRINT THIS
PAGE

Description

Premier Implant Cement is a eugenol-free, elastic polymer cement indicated for use with long-term temporaries and implant-retained crowns that may require future removal. The cement has a two-stage cure – an initial gel-phase in two minutes 30 seconds for removal of excess cement and a rigid final set. It is a neutral, light-pink gingival color. The kit includes a 5-ml dual-barrel, auto-mix syringe; 10 mixing tips; and one implant scaler instrument. Thirty consultants evaluated **Premier Implant Cement** in the cementation of 175 crowns. This product received a 91% clinical rating.

Suggested Retail Cost

- \$51.35/kit

More About: [Premier Dental Products Co,](#)
[Cement - Temporary, Cement,](#)

Editors' Note: *The original evaluation was published in THE DENTAL ADVISOR, Vol. 22, No. 1, January/February 2005.*

CE0145

[FAQ](#) [Copyright Policy](#) [Privacy Policy](#) [Terms & Conditions](#) [Contact](#)

Dental Consultants, Inc. 3110 West Liberty, Ann Arbor, Michigan 48103
(800) 347-1330 (734) 665-2020 (734) 665-1648 fax info@dentaladvisor.com
© 2007 Dental Consultants, Inc.